

Rat monoclonal [23G3] Anti-Mouse IgE epsilon chain (Biotin)
ab99585

製品の概要

| | |
|----------|---|
| 製品名 | Rat monoclonal [23G3] Anti-Mouse IgE epsilon chain (Biotin) |
| 由来種 | Rat |
| ターゲット生物種 | Mouse |
| 特異性 | Reacts with the epsilon heavy chain of BALB/c mouse IgE as demonstrated by ELISA. |
| アプリケーション | 適用あり: Flow Cyt, ELISA, ICC/IF |
| 標識 | Biotin |

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. Store In the Dark. |
| バッファー | pH: 7.4 Preservative: 0.1% Sodium azide Constituent: PBS |
| 精製度 | Affinity purified |
| ポリ/モノ | モノクローナル |
| クローン名 | 23G3 |
| アイソタイプ | IgG1 |
| 軽鎖の種類 | kappa |

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab99585の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------|
| | | |

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| Flow Cyt | | Use 1µg for 10 ⁶ cells. ab18403 - Rat monoclonal IgG1, is suitable for use as an isotype control with this antibody. |
| ELISA | | 1/5000 - 1/10000. |
| ICC/IF | | Use at an assay dependent concentration. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors