

Rat monoclonal [JC5-1] Anti-Mouse lambda light chain (Alkaline Phosphatase) ab99624

製品の概要

製品名	Rat monoclonal [JC5-1] Anti-Mouse lambda light chain (Alkaline Phosphatase)
由来種	Rat
ターゲット生物種	Mouse
特異性	ab99624 reacts with lambda 1 and lambda 2 light chains. It does not react with other immunoglobulin isotypes or kappa light chains.
アプリケーション	適用あり: ELISA, Flow Cyt, ICC/IF
標識	Alkaline Phosphatase

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 8.00 Preservative: 0.1% Sodium azide Constituents: 50% Glycerol (glycerin, glycerine), 0.0095% Magnesium chloride, 0.605% Tris
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	JC5-1
アイソタイプ	IgG2b
軽鎖の種類	kappa

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab99624の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		1/500 - 1/1000.
Flow Cyt		Use at an assay dependent dilution.
ICC/IF		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors