

Rat monoclonal [LO-MG1-2] Anti-Mouse IgG1 heavy chain (Biotin) ab11587

3 References

製品の概要

製品名	Rat monoclonal [LO-MG1-2] Anti-Mouse IgG1 heavy chain (Biotin)
由来種	Rat
ターゲット生物種	Mouse
特異性	ab11587 binds to the heavy chain of murine IgG1.
アプリケーション	適用あり: ELISA
免疫原	Full length protein corresponding to Mouse IgG1. Purified mouse IgG1 from BALB/c mice
標識	Biotin

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40 Preservative: 0.1% Sodium azide Constituent: PBS
精製度	Protein G purified
特記事項 (精製)	Purified IgG prepared by affinity chromatography.
ポリ/モノ	モノクローナル
クローン名	LO-MG1-2
ミエローマ	IR983F
アイソタイプ	IgG1
特記事項	This product should be stored undiluted. Storage in frost free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use.

アプリケーション

The Abpromise guarantee

Abpromise保証は、 次のテスト済みアプリケーションにおけるab11587の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		

追加情報

ELISA: Use at a concentration of 0.125 µg/ml, as detecting antibody.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors