

Rat monoclonal [LO-MM-3] Anti-Mouse IgM heavy chain ab11658

製品の概要

製品名	Rat monoclonal [LO-MM-3] Anti-Mouse IgM heavy chain
由来種	Rat
ターゲット生物種	Mouse
特異性	ab11658 recognises murine IgM (u heavy chain) and does not cross react with other murine immunoglobulin classes or subclasses. Avidity = $2 \times 10^8 M^{-1}$
アプリケーション	適用あり: ELISA
免疫原	Polyclonal mouse IgM.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Purified IgM
特記事項(精製)	Purified IgM prepared by affinity chromatography on mouse IgM/Sepharose 4B.
ポリ/モノ	モノクローナル
クローン名	LO-MM-3
ミエローマ	IR983F
アイソタイプ	IgM
特記事項	This product should be stored undiluted. Storage in frost free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use. This antibody recognises murine IgM (μ heavy chain) and does not cross react with other murine immunoglobulin classes or subclasses. Avidity = $2 \times 10^8 M^{-1}$

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab11658の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		

追加情報

ELISA: Use at a concentration of 5 µg/ml (coating).

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors