

Rat monoclonal [SB77e] Anti-Mouse IgG1 H&L (HRP) preadsorbed ab99603

5 References

製品の概要

製品名	Rat monoclonal [SB77e] Anti-Mouse IgG1 H&L (HRP) preadsorbed
由来種	Rat
ターゲット生物種	Mouse
特異性	ab99603 has been tested by ELISA against a panel of purified Mouse Immunoglobulins including IgM, IgG1, IgG2a, IgG2b, IgG3, and IgA. This product is positive for IgG1 and has minimal reactivity with the other Mouse Immunoglobulins.
アプリケーション	適用あり: ELISA
吸着処理血清	Goat, Human, Rabbit
標識	HRP

[more details](#)

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.40 Constituents: 50% Glycerol (glycerin, glycerine), 50% PBS
精製度	IgG fraction
ポリ/モノ	モノクローナル
クローン名	SB77e
アイソタイプ	IgG2b

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab99603の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		1/4000 - 1/8000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors