

Rat monoclonal [SB77e] Anti-Mouse IgG1 H&L (PE) preadsorbed ab99605

画像数 1

製品の概要

製品名	Rat monoclonal [SB77e] Anti-Mouse IgG1 H&L (PE) preadsorbed
由来種	Rat
ターゲット生物種	Mouse
特異性	ab99605 reacts with Mouse IgG1 as demonstrated by ELISA. Each batch of Rat monoclonal [SB77e] Anti-Mouse IgG1 H&L (Phycoerythrin) preadsorbed is tested to show minimal cross reactivity with IgM, IgA, and other subsets of IgG.
アプリケーション	適用あり: ELISA, Flow Cyt, ICC/IF
吸着処理血清	Goat, Human, Rabbit
標識	PE. Ex: 488nm, Em: 575nm

[more details](#)

製品の特性


製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	pH: 7.3 Preservative: 0.1% Sodium azide Constituents: Sucrose, PBS
精製度	Affinity purified
ポリ/モノ	モノクローナル
クローン名	SB77e
アイソタイプ	IgG2b

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab99605の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.
Flow Cyt		Use at an assay dependent dilution.
ICC/IF		Use 0.1µg for 10 ⁶ cells.

画像



Top: (Negative control) ab99605 being used against Mouse IgG1 in human blood at 0.1 µg/10⁶ cells.

Bottom: ab99605 staining human blood incubated with mouse anti-human CD3 (IgG1 isotype) at 0.1 µg/10⁶ cells.

Flow Cytometry - Rat monoclonal [SB77e] Anti-Mouse IgG1 H&L (Phycoerythrin) preadsorbed (ab99605)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors