

# Bradford Reagent ab119216

**31 References**    **画像数 4**

### 製品の概要


製品名	Bradford Reagent
特記事項	<p>This product is manufactured by Expedeon, an Abcam company, and was previously called Bradford Ultra - 500 ml. BFU05L is the same as the 500 ml size.</p> <p>Bradford Reagent ab119216 offers an improved formulation over classical Bradford: it has higher tolerance of detergent contamination of the protein samples.</p> <p>This Bradford Reagent shows excellent linearity for a defined range of protein concentrations and shows significantly less protein-to-protein variation than is observed with other Bradford-type Coomassie formulations.</p> <p>When Coomassie dye binds protein in an acidic medium, an immediate shift in absorption maximum occurs from 465 nm to 595 nm with a concomitant color change from brown to blue. Protein concentrations are estimated by reference to the absorbance obtained for a series of standard protein dilutions, which are assayed alongside samples with unknown protein concentrations.</p> <p>Manufactured under patent numbers: US 12/226,797 ; SG 201103010-3 ; SG 200807953-5 ; JP 2009-507182 ; IN 5722/CHENP/2008 ; GB 0608377.8 ; EP 10187154.9; EP 07733647.7 ; CN 200780021434.8 ; CA 2,649,835 ; AU 2007245443.</p> <p>Support on how to determine protein concentration can be found <a href="#">here</a>.</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at Room Temperature. Store at +4°C.
バッファー	Constituents: 3% Ethanol, 2.3% α-Cyclodextrin hydrate, 0.23% gamma-Cyclodextrin, 14% Phosphoric acid

### 画像


Protein quantitation of nuclear extracts with ab113474 (using ab119216 Optiblot Bradford reagent)



Copyright (c) 2015 Abcam plc


Bradford Reagent (ab119216)

Protein quantitation of nuclear extracts with ab113474 (using competitor's Bradford)




Copyright (c) 2015 Abcam plc

Bradford Reagent (ab119216)



Copyright (c) 2015 Abcam plc

Bradford Reagent (ab119216)



**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors