

### PC-12 nuclear extract lysate (forskolin stimulated) ab14885

#### 製品の概要

製品名	PC-12 nuclear extract lysate (forskolin stimulated)
特記事項	The protein content has been determined by a Bradford-based assay. The PC12 nuclear extracts have been collected after a 30-minute incubation with or without forskolin (20 $\mu$ M).

#### 製品の特性

Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.50 Constituents: 0.37% PBS, 0.0002% Sodium molybdate, 0.216% Beta glycerophosphate, 0.018% Sodium orthovanadate, 0.042% Sodium fluoride, 0.015% DTT, 0.48% HEPES, 0.003% EDTA, 20% Glycerol, 2.34% Sodium chloride
背景	With Protease inhibitor Generated by Greene and Tischler (1976) from a transplantable rat adrenal pheochromocytoma line, PC12 cells respond reversibly to nerve growth factor (NGF) and therefore have a neuronal phenotype, making the line a useful model system for neuronal differentiation. They produce catecholamines (dopamine, norepinephrine), neurotransmitters, NPP and are therefore used in studies of differentiation, neurobiochemistry and neurobiology.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors