


Regulation of Histone Modifications



Chromatin products from Abcam:

See more products at: www.abcam.com

Product	Clonality	Applications	Datasheet
			www.abcam.com/ab...
ChIP-on-chip Antibodies			
Histone H3 (acetyl K18)	P	ChIP, ChIP on chip, ICC/IF, IHC-P, WB	1191
Histone H3 (acetyl K27)	P	ChIP, ChIP on chip, IF, WB	4729
Histone H3 (di methyl K4)	P	ChIP, ChIP on chip, Dot, ICC/IF, IHC-P, WB	7766
Histone H3 (mono methyl K4)	P	ChIP, ChIP on chip, Flow Cyt., ICC/IF, IF, WB	8895
Histone H3 (mono methyl R2)	P	ChIP, ChIP on chip, WB	15584
Histone H3 (tri methyl K27)	M	ChIP, ChIP on chip, ICC/IF, WB	6002
[mAbcam 6002]			
Histone H3 (tri methyl K36)	P	ChIP, ChIP on chip, IHC-P, WB	9050
Histone H3 (tri methyl K4)	P	ChIP, ChIP on chip, Flow Cyt., ICC, ICC/IF, IF, WB	8580
Histone H3 (tri methyl K79)	P	ChIP, ChIP on chip, WB	2621
Histone H3 (tri methyl K9)	P	ChIP, ChIP on chip, ICC/IF, WB	8898
Histone H3	P	ChIP, ChIP on chip, ICC, ICC/IF, IF, IP, WB	1791
Histone H4 (acetyl K5)	P	ChIP, ChIP on chip, ELISA, IF, IP, WB	1758
Histone ChIP grade Antibodies			
Histone H1.2	P	ChIP, ICC/IF, WB	4086
Histone H2A (acetyl K5)	P	ChIP, ELISA, ICC/IF, IF, IP, WB	1764
Histone H2A (phospho S129)	P	ChIP, WB	15083
Histone H2A	P	ChIP, WB	15653
Histone H2A	P	ChIP, ICC/IF, WB	18255
Histone H2A,Z	P	ChIP, IF, WB	18263
Histone H2A,Z	P	ChIP, ICC/IF, WB	4174
Histone H2B (acetyl K5) [EP857Y]	M	ChIP, ICC/IF, IHC-P, IP, WB	40886
Histone H2B (acetyl)	P	ChIP, ELISA, ICC/IF, IF, IP, WB	1759
Histone H2B	P	ChIP, ICC/IF, WB	1790
Histone H3 (acetyl K9)	P	ChIP, WB	10812
Histone H3 (acetyl K9)	P	ChIP, Dot, ICC/IF, WB	4441
Histone H3 (acetyl K8)	P	ChIP, ICC/IF, IF, WB	12178
Histone H3 (acetyl K9) [AH3-120]	M	ChIP, ELISA, ICC, ICC/IF, IF, WB	12179
Histone H3 (asymmetric di methyl R17)	P	ChIP, Dot, IF, IHC-P, IP, WB	8284
Histone H3 (asymmetric di methyl R2)	P	ChIP, IF, WB	8046
Histone H3 (asymmetric di methyl R26)	P	ChIP, IF, WB	8047
Histone H3 (citulline 2 + 8 + 17) [CH3]	P	ChIP, IF, WB	5103
Histone H3 (di methyl K27)	P	ChIP, WB	24684
Histone H3 (di methyl K4)	P	ChIP, WB	11946
Histone H3 (di methyl K4) [Y47]	M	ChIP, ICC/IF, IHC-P, WB	32356
Histone H3 (di methyl K79)	P	ChIP, WB	3594
Histone H3 (di methyl K9)	M	ChIP, ICC/IF, WB	1220
[mAbcam 1220]			
Histone H3 (di+tri methyl K4)	M	ChIP, ICC, IF, WB	6000
[mAbcam 6000]			
Histone H3 (methyl K79)	P	ChIP, WB	28940

Product	Clonality	Applications	Datasheet
			www.abcam.com/ab...
Histone ChIP grade Antibodies			
Histone H3 (mono methyl K27, tri methyl K27 + K4) [8H2]	M	ChIP, IF, IP, WB	14222
Histone H3 (mono methyl K36)	P	ChIP, ICC, IF, WB	9048
Histone H3 (mono methyl K79)	P	ChIP, ICC, IF, WB	2886
Histone H3 (mono methyl K9)	P	ChIP, WB	8896
Histone H3 (mono methyl K9)	P	ChIP, Flow Cyt., IF, WB	9045
[mAbcam 14955]			
Histone H3 (phospho S10)	M	ChIP, ICC/IF, WB	14955
Histone H3 (tri methyl K4)	M	ChIP, ICC, IF, WB	1012
[mAbcam 1012]			
Histone H3 (tri methyl K4)	M	ChIP, WB	12209
[mAbcam 12209]			
Histone H3	P	ChIP, WB	12079
Histone H3	P	ChIP, ELISA, IF, IP, WB	7834
Histone H3 [mAbcam 10799]	M	ChIP, WB	10799
Histone H4 (acetyl K12)	P	ChIP, ELISA, ICC/IF, IHC-P, IP, WB	1761
Histone H4 (acetyl K16)	P	ChIP, ELISA, ICC/IF, IP, WB	1762
Histone H4 (acetyl K8)	P	ChIP, ELISA, IF, IP, WB	1760
Histone H4 (acetyl K91)	P	ChIP, WB	4627
Histone H4 (mono methyl K20)	P	ChIP, Flow Cyt., IF, WB	9051
Histone H4 (symmetric di methyl R3)	P	ChIP, ICC/IF, WB	5823
Histone H4 (tri methyl K20)	P	ChIP, Flow Cyt., IF, WB	9053
Histone H4 [mAbcam31827]	M	ChIP, WB	31827
Novel Antibodies			
Histone H2A (asymmetric di methyl R3)	P	Fast track	21574
Histone H2A (Hydroxy P26)	P	Fast track	30819
Histone H2A (mono methyl K125)	P	Fast track	23633
Histone H2A (mono methyl K127)	P	Fast track	12950
Histone H2A (mono methyl R17)	P	Fast track	28002
Histone H2A (tri methyl K125)	P	Fast track	23958
Histone H2A (tri methyl K127)	P	Fast track	23957
Histone H2B (di methyl K43)	P	Fast track	4629
Histone H2B (Hydroxy P10)	P	Fast track	30865
Histone H3 (Hydroxy P121)	P	Fast track	26103
Histone H3 (Hydroxy P38)	P	Fast track	19849
Histone H4 (Hydroxy P32)	P	Fast track	26104
Histone H4 (tri methyl K59)	P	Fast track	21084
ChIP related products			
Rabbit Control IgG	P	ChIP, ELISA, WB	46540
ss Deoxyribonucleic acid from herring sperm	M	BL, ChIP, ISH	46666

Active

Intergenic

↑

Gene

Intergenic

View all our Histone Modification antibodies at www.abcam.com