






DNA Repair Pathways



www.abcam.com/nuclearsignal

Abcam in the USA
Abcam Inc.
1 Kendall Square, Ste 341
Cambridge,
MA 02139-1517 USA

Tel: 1-617-225-2272
Toll free: 1-888-77-ABCAM
Toll free Fax: 1-866-739-9884

Abcam in Europe
Abcam plc,
330 Cambridge Science Park
Milton Road,
Cambridge,
CB4 0FL UK

Tel: +44 1223 696000
Fax: +44 1223 771600

Abcam in Japan
Abcam K.K.
1-16-8 Nihonbashi Kakigaracho
Chuo-ku,
Tokyo 103-0014
Japan

Tel: +81-(0)3-6231-0940
Fax: +81-(0)3-6231-0941

Featured DNA Repair antibodies

| | | | | | | | |
|--|---------------------------|--|------------------------------|---|----------------------------|---|---|
| | ab29 PCNA antibody | | ab32138 PARP antibody | | ab2888 ATM antibody | | ab2893 gamma H2A.X (phospho S139) antibody |
| Formaldehyde fixed mouse tissue sections stained with ab29 at 1/6000 dilution. | | Paraffin-embedded sections of human brain stained with ab32138 at 1/25 dilution. | | HeLa cells stained with ab2888 at 1/200 dilution. | | Pachytene spermatocytes stained with ab2893 at 1/200. | |