

Product datasheet

2X LDL/VLDL Precipitation Buffer ab105138

1 References

製品の概要

製品名	2X LDL/VLDL Precipitation Buffer
アッセイタイプ	Quantitative
製品の概要	2X LDL/VLDL Precipitation Buffer for ab65390 HDL and LDL/VLDL Cholesterol Assay Kit .

製品の特性

保存方法 Store at -20°C.

内容	ラベル	10 ml
2X LDL/VLDL Precipitation Buffer	NM	1 x 10ml

関連性 HDL (high-density-lipoprotein), LDL (low-density lipoprotein) and VLDL (very low-density lipoprotein) are lipoproteins that act as carrier proteins for cholesterol. The regulation of HDL-cholesterol and LDL-cholesterol plays a central role in various disease development. It is well known that low levels of HDL and high level of LDL are associated with an increased risk of cardiovascular events.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors