

Product datasheet

Anti-ABCA5 antibody ab82582

製品の概要

製品名 Anti-ABCA5 antibody  
 製品の詳細 Goat polyclonal to ABCA5

**i** This product is a **fast track antibody**. It has been affinity purified and shows high titre values against the immunizing peptide by ELISA.

[Read the terms of use »](#)

由来種 Goat

種交差性

交差が予測される動物種: Human 

免疫原

Synthetic peptide:  
 DYSSETSEDDDSLKC  
 , corresponding to internal sequence amino acids 1361 - 1375 of human ABCA5 (NP\_758424.1)

 [Run BLAST with](#)

 [Run BLAST with](#)

製品の特性

製品の状態 Liquid  
 保存方法 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.  
 バッファー pH: 7.30  
 Preservative: 0.02% Sodium azide  
 Constituents: 0.5% Tris buffered saline, 0.5% BSA  
 精製度 Immunogen affinity purified  
 特記事項 (精製) ab82582 is purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.  
 ポリ/モノ ポリクローナル  
 アイソタイプ IgG

アプリケーション

追加情報 Peptide ELISA: antibody detection limit dilution 1:1000.  
 WB: This antibody did not give bands in initial tests on lysate/cell line. We would welcome your

feedback if you have tested this antibody in more suitable positive controls.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

---

機能	May play a role in the processing of autolysosomes.
組織特異性	Ubiquitously expressed. Highly expressed in testis, skeletal muscle, kidney, liver and placenta.
配列類似性	Belongs to the ABC transporter superfamily. ABCA family. Contains 2 ABC transporter domains.
発生段階	Expressed in fetal liver, kidney and brain.
翻訳後修飾	N-glycosylated.
細胞内局在	Lysosome membrane. Late endosome membrane. Golgi apparatus membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors