

Anti-AGTRAP antibody ab85175

面/色粒?

医薬用外毒物

製品の概要

製品名	Anti-AGTRAP antibody
製品の詳細	Rabbit polyclonal to AGTRAP
由来種	Rabbit
アプリケーション	適用あり: IHC-P, WB
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to Human AGTRAP (C terminal). Database link: NM_020350
ポジティブ・コントロール	WB: HeLa, HepG2, SK0V3, MCF-7 and MDA-MD-231 cell lysates. IHC-P: Human kidney tissue.

製品の特性

製品の状態	Lyophilised:Reconstitute in 200ul sterile H2O.
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituent: 1% BSA
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab85175** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/200.

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB 1/200 - 1/1000. Predicted molecular weight: 17 kDa.

ターゲット情報

画像



Western blot - Anti-AGTRAP antibody (ab85175)

All lanes : Anti-AGTRAP antibody (ab85175) at 1/500 dilution

Lane 1 : HeLa cell lysate


Lane 2 : HepG2 cell lysate

Lane 3 : SKOV3 cell lysate

Lane 4 : MCF-7 cell lysate

Lane 5 : MDA-MD-231 cell lysate

Predicted band size: 17 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-AGTRAP antibody (ab85175)

ab85175, at 1/100 dilution, staining AGTRAP in formalin-fixed, paraffin-embedded human kidney by Immunohistochemistry. Note cytoplasmic staining.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors