

## Product datasheet

# Anti-ANXA9 antibody ab37501

### 画像数 2

#### 製品の概要

製品名	Anti-ANXA9 antibody
製品の詳細	Chicken polyclonal to ANXA9
由来種	Chicken
アプリケーション	<b>適用あり:</b> ELISA, WB
種交差性	<b>交差種:</b> Human
免疫原	Full length native protein (purified), corresponding to amino acids 1-338 of Human ANXA9

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgY

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab37501** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use a concentration of 1 µg/ml.
WB		1/500. Predicted molecular weight: 38 kDa.

#### ターゲット情報


<b>機能</b>	Low affinity receptor for acetylcholine known to be targeted by disease-causing pemphigus vulgaris antibodies in keratinocytes.
<b>組織特異性</b>	Expressed in the stratified squamous skin epithelium, but not in epithelia of other types (at protein level).
<b>配列類似性</b>	Belongs to the annexin family. Contains 4 annexin repeats.

## 画像



ELISA - Anti-ANXA9 antibody (ab37501)

Direct ELISA test. Native Human Annexin A9 was test antigen. ab37501 (chicken IgY anti-annexin A9) was primary antibody and Goat anti-IgY-HRP was secondary antibody. A fixed amount of antibody (1 µg/ml) and serial dilutions of antigen were used.



Western blot - Anti-ANXA9 antibody (ab37501)

Anti-ANXA9 antibody (ab37501) at 1/2000 dilution + E. coli-derived fusion protein

### Secondary

Goat anti-IgY-HRP at 1/1000 dilution

**Predicted band size:** 38 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors