


Product datasheet

Anti-APBB3 antibody ab96814

画像数 1

製品の概要

製品名	Anti-APBB3 antibody
製品の詳細	Rabbit polyclonal to APBB3
由来種	Rabbit
特異性	ab96814 recognizes APBB3 Isoform I.
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Cow 
免疫原	Recombinant protein fragment corresponding to a region within amino acids 240 and 483 of Human APBB3 (NP_573420; UniProt ID: O95704 isoform I)
ポジティブ・コントロール	A431 whole cell lysate; 293T, H1299, MOLT4 and Raji cell lysates

法規制情報

医薬用外毒物

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 1.21% Tris, 0.75% Glycine, 20% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab96814** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 52 kDa.

ターゲット情報

機能	May modulate the internalization of beta-amyloid precursor protein.
組織特異性	Expressed in various tissues.
配列類似性	Contains 2 PID domains. Contains 1 WW domain.

画像



Anti-APBB3 antibody (ab96814) at 1/1000 dilution + A431 whole cell lysate at 30 µg

Predicted band size: 52 kDa

10% SDS Page

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors