

Product datasheet

CXCR6 overexpression 293T lysate (whole cell) ab94068

画像数 2

製品の概要

製品名	CXCR6 overexpression 293T lysate (whole cell)
製品の詳細	CXCR6 overexpression lysate
特記事項	ab94068 is a 293T cell transfected lysate in which Human Bonzo has been transiently over-expressed using a pCMV-Bonzo plasmid. The lysate is provided in 1X Sample Buffer. Note: For more detailed how the transfected lysate was prepared view preparation notes
アプリケーション	適用あり: WB

法規制情報

医薬用外劇物

製品の特性


製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituent: 100% 1x Sample Buffer
精製度	Whole Cell Lysate
背景	Function: Receptor for the C-X-C chemokine CXCL16. Used as a coreceptor by SIMs and by strains of HIV-2 and m-tropic HIV-1. Tissue specificity: Expressed in lymphoid tissues and activated T cells. Similarity: Belongs to the G-protein coupled receptor 1 family.

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab94068** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.



ab94068 at 15µg/lane on an SDS-PAGE gel.



All lanes : Anti-Bonzo antibody (ab89579) at 1/500 dilution

Lane 1 : CXCR6 overexpression 293T lysate (whole cell) (ab94068) at 25 µg

Lane 2 : 293T non-transfected lysate

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors