

Product datasheet

Anti-C18orf8 antibody ab78596

画像数 2

製品の概要

製品名	Anti-C18orf8 antibody
製品の詳細	Goat polyclonal to C18orf8
由来種	Goat
アプリケーション	適用あり: IHC-P, WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Dog 
免疫原	Synthetic peptide: C-DKLNHEYKKYLDAEQ , corresponding to internal sequence amino acids 414-428 of Human C18orf8 (NP_037458.3) Run BLAST with Run BLAST with
ポジティブ・コントロール	Human bone marrow, spleen or tonsil lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 0.5% BSA, 0.5% Tris buffered saline
精製度	Immunogen affinity purified
特記事項 (精製)	Ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab78596** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 5 - 10 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		Use a concentration of 1 - 3 µg/ml. Detects a band of approximately 80 kDa (predicted molecular weight: 75 kDa).

ターゲット情報

配列類似性


Belongs to the Mic1 family.
Contains 1 Mic1 domain.

画像



ab87596 (5µg/ml) staining of C18orf8 in paraffin embedded Human Small Intestine following steamed antigen retrieval with citrate buffer pH 6 and AP-staining shows textured cytoplasm staining in the enterocytes.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-C18orf8 antibody (ab78596)



Anti-C18orf8 antibody (ab78596) at 1/1 dilution + Human breast tissue lysate in RIPA buffer at 35 µg

Developed using the ECL technique.

Predicted band size: 75 kDa

Observed band size: 80 kDa

Western blot - Anti-C18orf8 antibody (ab78596)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors