

## Product datasheet

# Anti-C3 antibody ab97462

★★★★★ 1 Abreviews 5 References 画像数 7

### 製品の概要

製品名	Anti-C3 antibody
製品の詳細	Rabbit polyclonal to C3
由来種	Rabbit
アプリケーション	適用あり: IHC-P, IP, WB, ICC/IF
種交差性	交差種: Human 交差が予測される動物種: Mouse
免疫原	Recombinant fragment, corresponding to a region within amino acids 1498 - 1625 of Human C3.
ポジティブ・コントロール	WB: HepG2 and Huh7 whole cell lysate. ICC/IF: HeLa cells. IHC: H1299 xenograft tissue; mouse brain tissue. IP: HepG2 whole cell extract.

### 法規制情報

**医薬用外毒物**

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 40% Glycerol, 1X PBS, pH 7.0
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab97462** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P	★★★★★	1/100 - 1/1000.
IP		1/100 - 1/500.
WB		1/500 - 1/3000. Predicted molecular weight: 187 kDa.
ICC/IF		1/100 - 1/1000.

## ターゲット情報


<b>機能</b>	<p>C3 plays a central role in the activation of the complement system. Its processing by C3 convertase is the central reaction in both classical and alternative complement pathways. After activation C3b can bind covalently, via its reactive thioester, to cell surface carbohydrates or immune aggregates.</p> <p>Derived from proteolytic degradation of complement C3, C3a anaphylatoxin is a mediator of local inflammatory process. It induces the contraction of smooth muscle, increases vascular permeability and causes histamine release from mast cells and basophilic leukocytes.</p>
<b>組織特異性</b>	Plasma.
<b>関連疾患</b>	<p>Defects in C3 are the cause of complement component 3 deficiency (C3D) [MIM:120700]. A rare defect of the complement classical pathway. Patients develop recurrent, severe, pyogenic infections because of ineffective opsonization of pathogens. Some patients may also develop autoimmune disorders, such as arthralgia and vasculitic rashes, lupus-like syndrome and membranoproliferative glomerulonephritis.</p> <p>Genetic variation in C3 is associated with susceptibility to age-related macular degeneration type 9 (ARMD9) [MIM:611378]. ARMD is a multifactorial eye disease and the most common cause of irreversible vision loss in the developed world. In most patients, the disease is manifest as ophthalmoscopically visible yellowish accumulations of protein and lipid that lie beneath the retinal pigment epithelium and within an elastin-containing structure known as Bruch membrane.</p> <p>Defects in C3 are a cause of susceptibility to hemolytic uremic syndrome atypical type 5 (AHUS5) [MIM:612925]. An atypical form of hemolytic uremic syndrome. It is a complex genetic disease characterized by microangiopathic hemolytic anemia, thrombocytopenia, renal failure and absence of episodes of enterocolitis and diarrhea. In contrast to typical hemolytic uremic syndrome, atypical forms have a poorer prognosis, with higher death rates and frequent progression to end-stage renal disease. Note=Susceptibility to the development of atypical hemolytic uremic syndrome can be conferred by mutations in various components of or regulatory factors in the complement cascade system. Other genes may play a role in modifying the phenotype.</p>
<b>配列類似性</b>	<p>Contains 1 anaphylatoxin-like domain.</p> <p>Contains 1 NTR domain.</p>
<b>翻訳後修飾</b>	<p>C3b is rapidly split in two positions by factor I and a cofactor to form iC3b (inactivated C3b) and C3f which is released. Then iC3b is slowly cleaved (possibly by factor I) to form C3c (beta chain + alpha' chain fragment 1 + alpha' chain fragment 2), C3dg and C3f. Other proteases produce other fragments such as C3d or C3g.</p> <p>Phosphorylation sites are present in the extracellular medium.</p>
<b>細胞内局在</b>	Secreted.



Immunocytochemistry/ Immunofluorescence - Anti-C3 antibody (ab97462)

HeLa (human epithelial cell line from cervix adenocarcinoma) cells labeling C3 with ab97462 at 1/200 dilution (green) in ICC/IF. Cells were fixed in 4% paraformaldehyde at room temperature for 15 minutes.

Nuclei were stained using Hoechst 33342 (blue).



Western blot - Anti-C3 antibody (ab97462)

**All lanes** : Anti-C3 antibody (ab97462) at 1/1000 dilution

**Lane 1** : Huh7 whole cell lysate

**Lane 2** : HepG2 (human liver hepatocellular carcinoma cell line) whole cell lysate


Lysates/proteins at 30 µg per lane.

**Secondary**

**All lanes** : HRP-conjugated anti-rabbit IgG


**Predicted band size:** 187 kDa

5% SDS-PAGE.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-C3 antibody (ab97462)

Paraffin-embedded mouse brain tissue stained for C3 using ab97462 at 1/500 dilution in immunohistochemical analysis.




Immunoprecipitation - Anti-C3 antibody (ab97462)

C3 was immunoprecipitated from HepG2 (human liver hepatocellular carcinoma cell line) whole cell extracts using 5  $\mu$ g of ab97462. Western blot was performed from the immunoprecipitate using ab97462.

**Lane 1:** HepG2 whole cell extract.

**Lane 2:** Control IgG instead of ab97462 in HepG2 whole cell extract.

**Lane 3:** ab97462 IP in HepG2 whole cell extract.




Western blot - Anti-C3 antibody (ab97462)

Anti-C3 antibody (ab97462) at 1/1500 dilution + HepG2 whole cell lysate at 30  $\mu$ g


**Predicted band size:** 187 kDa

7.5% SDS-PAGE.



Immunocytochemistry/ Immunofluorescence - Anti-C3 antibody (ab97462)

ab97462 at 1/100 dilution staining C3 in HeLa cells by Immunofluorescence, Paraformaldehyde fixed. Lower image shows cells co-stained with Hoechst 33342.



ab97462 at 1/100 dilution staining C3 in H1299 xenograft by Immunohistochemistry, Paraffin-embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-C3 antibody (ab97462)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors