

Product datasheet

Anti-C3a/C3a(desArg)/C3 antibody [D17/1] (Biotin) ab121050

画像数 1

製品の概要

製品名	Anti-C3a/C3a(desArg)/C3 antibody [D17/1] (Biotin)
製品の詳細	Mouse monoclonal [D17/1] to C3a/C3a(desArg)/C3 (Biotin)
由来種	Mouse
標識	Biotin
アプリケーション	適用あり: Sandwich ELISA, IHC-P
種交差性	交差種: Human
免疫原	Human C3a/C3a des Arg.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40 Preservative: 0.1% Sodium azide Constituents: 99% Phosphate Buffer, 0.82% Sodium chloride
精製度	Protein G purified
特記事項(精製)	Protein-A/G purified.
ポリ/モノ	モノクローナル
クローン名	D17/1
アイソタイプ	IgG1
軽鎖の種類	kappa


アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab121050** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Sandwich ELISA		1/40000. Can be paired for Sandwich ELISA with Mouse monoclonal [K13/16] to C3a / C3a des Arg (ab36385) , as detection antibody. The dilution guideline for ELISA is based on use as detection antibody for C3a / C3a des Arg coated at 1 µg/ml.
IHC-P		Use at an assay dependent concentration.

画像



Immunohistochemistry of human appendix section staining C3a/C3a(desArg)/C3 with ab121050 (1:30)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-C3a/C3a(desArg)/C3 antibody [D17/1] (Biotin) (ab121050)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors