

## Product datasheet

# Anti-C5a / C5a des Arg antibody [2952] ab11878

★★★★☆ 1 Abreviews 3 References 画像数 1

### 製品の概要

製品名	Anti-C5a / C5a des Arg antibody [2952]
製品の詳細	Mouse monoclonal [2952] to C5a / C5a des Arg
由来種	Mouse
特異性	This antibody reacts with a neo-epitope (des-Arg) on C5a that is formed upon cleavage of C5 into C5a and C5b. It will detect C5a and C5a des-Arg (Not C5).
アプリケーション	<b>適用あり:</b> IHC-Fr, Functional Studies, IHC-P, ELISA, WB
種交差性	<b>交差種:</b> Human
免疫原	ISHKDMQLG (C5a amino acid 65-73)
	<a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
ポジティブ・コントロール	Human Granulocytes

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: PBS, 0.1% BSA
精製度	Protein G purified
特記事項(精製)	0.2 µm filtered
ポリ/モノ	モノクローナル
クローン名	2952
アイソタイプ	IgG1

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab11878** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-Fr		Use at an assay dependent concentration.
Functional Studies		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
ELISA	★★★★☆	Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

## ターゲット情報

**関連性** C5 is involved in the activation of the lytic pathway within the complement system, which is an important factor in innate immunity.

## 画像



ELISA with [ab11876](#) as capture antibody at 1/100, ab11878 as detection antibody at 1/250, diluted in carbonate buffer or PBS. Dilutions of commercial C5a were used as the standard.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors