

Product datasheet

Anti-C6orf62 antibody ab107371

画像数 3

製品の概要

製品名	Anti-C6orf62 antibody
製品の詳細	Rabbit polyclonal to C6orf62
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Mouse, Hamster, Human
免疫原	Synthetic peptide conjugated to KLH, between amino acids 118-147 from the central region of Human C6orf62 (NP_112201.1).
ポジティブ・コントロール	CHO cell line lysate, Mouse kidney tissue lysate, Human brain tissue

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
特記事項(精製)	ab107371 is purified through a protein A column, followed by peptide affinity purification.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab107371** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/100 - 1/500. Predicted molecular weight: 27 kDa.

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.


## ターゲット情報

## 画像




Anti-C6orf62 antibody (ab107371) at 1/100 dilution + CHO cell line lysate at 35 µg

**Predicted band size:** 27 kDa



Anti-C6orf62 antibody (ab107371) at 1/100 dilution + Mouse kidney tissue lysate at 35 µg

**Predicted band size:** 27 kDa



Immunohistochemistry analysis in formalin fixed and paraffin embedded Human brain tissue followed by peroxidase conjugation of the secondary antibody and DAB staining. ab107371 at a dilution of 1/50.

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors