

Product datasheet

Anti-CCDC153 antibody ab111214

画像数 1

製品の概要

製品名	Anti-CCDC153 antibody
製品の詳細	Rabbit polyclonal to CCDC153
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a regionn within N terminal amino acids 39-68 of Human CCDC153 (NP_001138490.1)
ポジティブ・コントロール	HepG2 cell line lysates
特記事項	Prevent freeze-thaw cycles.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
特記事項(精製)	This antibody is purified through a protein A column, followed by peptide affinity purification.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab111214** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/100 - 1/500. Predicted molecular weight: 24 kDa.

## ターゲット情報

---


### 配列類似性

Belongs to the UPF0610 family.

---

## 画像

---



Western blot - Anti- CCDC153 antibody (ab111214)

Anti-CCDC153 antibody (ab111214) at 1/100 dilution + HepG2 cell line lysates at 35 µg

**Predicted band size:** 24 kDa

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors