

## Product datasheet

# Anti-CCDC83 antibody ab97917

画像数 2

### 製品の概要

製品名	Anti-CCDC83 antibody
製品の詳細	Rabbit polyclonal to CCDC83
由来種	Rabbit
アプリケーション	適用あり: WB, ICC/IF
種交差性	交差種: Human 交差が予測される動物種: Cow
免疫原	Recombinant protein fragment, corresponding to a region within amino acids 17-276 of Human CCDC83 (AAH40280).
ポジティブ・コントロール	293T, A431, H1299, HeLa, HepG2, MOLT4 and Raji cell lysates for WB; HeLa cells for IF/ICC

### 法規制情報

**医薬用外毒物**

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 0.75% Glycine, 10% Glycerol, 1.21% Tris
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab97917** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 49 kDa.
ICC/IF		1/100 - 1/200.

## ターゲット情報

**関連性** There are three named isoforms.


## 画像



Anti-CCDC83 antibody (ab97917) at 1/1000 dilution + HepG2 whole cell lysate at 30 µg

**Predicted band size:** 49 kDa

10% SDS PAGE



Immunofluorescence analysis of paraformaldehyde-fixed HeLa cells, using ab97917 at 1/200 dilution. Lower image is merged with DNA probe.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors