




## Product datasheet

# Anti-CD134 / OX40L receptor antibody ab111475

画像数 1

### 製品の概要

製品名	Anti-CD134 / OX40L receptor antibody
製品の詳細	Rabbit polyclonal to CD134 / OX40L receptor
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Pig 
免疫原	Synthetic peptide within Human CD134/ OX40L receptor aa 151-200 (internal sequence). The exact sequence is proprietary. Sequence: GKHTLQPASNSSDAICEDRDPPATQPQETQGPPARPITVQ PTEAWPRTSQ  Database link: <a href="#">NP_003318</a>   <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a>
ポジティブ・コントロール	WB: 721_B cell lysate.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab111475** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.2 - 1 µg/ml. Predicted molecular weight: 29 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## ターゲット情報

機能	Receptor for TNFSF4/OX40L/GP34.
配列類似性	Contains 4 TNFR-Cys repeats.
細胞内局在	Membrane.

## 画像



Anti-CD134 / OX40L receptor antibody (ab111475) at 1 µg/ml +  
721\_B cell lysate at 10 µg

**Predicted band size:** 29 kDa

Gel concentration: 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors