

Product datasheet

Anti-CD2 antibody ab37212

1 References [画像数 2](#)

製品の概要

製品名	Anti-CD2 antibody
製品の詳細	Rabbit polyclonal to CD2
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Recombinant fragment corresponding to CD2.
ポジティブ・コントロール	Tonsil.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.40 Preservative: 0.1% Sodium azide Constituents: 1% BSA, 0.0268% PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab37212** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration. Detects a band of approximately 55 kDa (predicted molecular weight: 39 kDa).

アプリケーション	Abreviews	特記事項
IHC-P		1/50. Incubate for 10 min at RT. Staining of formalin-fixed tissues requires boiling tissue sections in 1mM EDTA, pH 8.0 for 10 min followed by cooling at RT for 20 min.

ターゲット情報


機能	CD2 interacts with lymphocyte function-associated antigen (LFA-3) and CD48/BCM1 to mediate adhesion between T-cells and other cell types. CD2 is implicated in the triggering of T-cells, the cytoplasmic domain is implicated in the signaling function.
配列類似性	Contains 1 Ig-like C2-type (immunoglobulin-like) domain. Contains 1 Ig-like V-type (immunoglobulin-like) domain.
細胞内局在	Membrane.

画像



Ab37212 staining human CD2 in tonsil by immunohistochemistry, Formalin fixed Paraffin embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD2 antibody (ab37212)



Anti-CD2 antibody (ab37212) at 1/25 dilution
+ Jurkat cell lysate

Predicted band size: 39 kDa

Observed band size: 55 kDa

Western blot - Anti-CD2 antibody (ab37212)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors