

Product datasheet

Anti-CD229 antibody [HLy9.25], prediluted (Phycoerythrin)
ab103171

1 References [画像数 1](#)

製品の概要

製品名	Anti-CD229 antibody [HLy9.25], prediluted (Phycoerythrin)
製品の詳細	Mouse monoclonal [HLy9.25] to CD229, prediluted (Phycoerythrin)
由来種	Mouse
標識	Phycoerythrin. Ex: 488nm, Em: 575nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Human
免疫原	CD299-transfected 300.19 pre-B cell line (Tissue/ cell preparation)
特記事項	Clone HLy9.25 is also known as HLy9.1.25.

製品の特性

製品の状態	Prediluted
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.097% Sodium azide Constituents: PBS, 0.2% BSA
精製度	Size exclusion
特記事項 (精製)	Conjugates are made from highly purified antibodies (>95% by SDS-PAGE).
ポリ/モノ	モノクローナル
クローン名	HLy9.25
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab103171** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 20µl for 10 ⁶ cells. or 100 µl of whole blood

ターゲット情報

機能	May participate in adhesion reactions between T lymphocytes and accessory cells by homophilic interaction.
配列類似性	Contains 2 Ig-like C2-type (immunoglobulin-like) domains. Contains 2 Ig-like V-type (immunoglobulin-like) domains.
細胞内局在	Membrane.

画像



Flow Cytometry analysis of labeling CD229 with Anti-CD229 antibody [HLy9.25], prediluted (Phycoerythrin) (ab103171).

Flow Cytometry - Anti-CD229 antibody [HLy9.25], prediluted (Phycoerythrin) (ab103171)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors