

Anti-CD23 antibody [B-G6] (Phycoerythrin) ab27304

画像 1

医薬用外毒物

製品の概要

製品名	Anti-CD23 antibody [B-G6] (Phycoerythrin)
製品の詳細	Mouse monoclonal [B-G6] to CD23 (Phycoerythrin)
由来種	Mouse
標識	Phycoerythrin. Ex: 488nm, Em: 575nm
特異性	This antibody is specific for CD23
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Human
免疫原	Daudi cell line
特記事項	Abcam is committed to meeting high standards of ethical manufacturing and has decided to discontinue this product by June 2019 as it has been generated by the ascites method. We are sorry for any inconvenience this may cause.

製品の特性

製品の状態	Lyophilised:Reconstitute with 1 ml deionised water
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.1% Sodium azide Constituents: 1% PBS, 1% BSA
精製度	Ion Exchange Chromatography
ポリ/モノ	モノクローナル
クローン名	B-G6
ミエローマ	x63-Ag8.653
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab27304** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

Flow Cyt

追加情報

Flow Cyt: Use 10µl for 10⁶ cells or 100 µl of whole blood.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

機能

This receptor has essential roles in the regulation of IgE production and in the differentiation of B-cells (it is a B-cell-specific antigen).

配列類似性

Contains 1 C-type lectin domain.


翻訳後修飾

N- and O-glycosylated.

細胞内局在

Cell membrane. Secreted. Also exists as a soluble excreted form.

画像



A typical staining pattern with ab27304 (10 µl per 10⁶ cells or 100 µl of whole blood) of lymphocytes.

Flow Cytometry - Anti-CD23 antibody [B-G6]
(Phycoerythrin) (ab27304)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors