

Product datasheet

Anti-CHD3 antibody ab84528

★★★★☆ 1 Abreviews 2 References 画像数 3

製品の概要

製品名	Anti-CHD3 antibody
製品の詳細	Rabbit polyclonal to CHD3
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Mouse, Human 交差が予測される動物種: Rat, Rabbit, Horse, Guinea pig, Dog, Pig, Chimpanzee, Rhesus monkey, Gorilla, Chinese hamster, Orangutan ▲
免疫原	Synthetic peptide corresponding to a region between residues 1950 and 2000 of Human CHD3 (NP_001005273.1).
ポジティブ・コントロール	Human Anaplastic Thyroid Carcinoma, Ewing Sarcoma, Non-Small Cell Lung Cancer, Ovarian Carcinoma, Pancreatic Islet Cell Tumor, Prostate Carcinoma, Skin Basal Cell Carcinoma, Skin Squamous Cell Carcinoma and Small Cell Lung Cancer tissue. Mouse Teratoma and Renal Cell Carcinoma tissue.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 0.1% BSA, Tris buffered saline
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab84528** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		1/100 - 1/500. Antigen retrieval is not essential but may optimise staining. The use of citrate buffer will enhance staining.

ターゲット情報


機能	Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones. Required for anchoring centrosomal pericentrin in both interphase and mitosis, for spindle organization and centrosome integrity.
組織特異性	Widely expressed.
配列類似性	<p>Belongs to the SNF2/RAD54 helicase family.</p> <p>Contains 2 chromo domains.</p> <p>Contains 1 helicase ATP-binding domain.</p> <p>Contains 1 helicase C-terminal domain.</p> <p>Contains 2 PHD-type zinc fingers.</p>
細胞内局在	Nucleus. Cytoplasm > cytoskeleton > centrosome. Associates with centrosomes in interphase and mitosis.

画像




ab84528, at 1/250 dilution, staining CHD3 in formalin-fixed, paraffin-embedded human prostate carcinoma by Immunohistochemistry using DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CHD3 antibody (ab84528)



ab84528, at 1/250 dilution, staining CHD3 in formalin-fixed, paraffin-embedded human basal cell carcinoma by Immunohistochemistry using DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CHD3 antibody (ab84528)



ab84528, at 1/250 dilution, staining CHD3 in formalin-fixed, paraffin-embedded mouse renal cell carcinoma by Immunohistochemistry using DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CHD3 antibody (ab84528)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors