

Product datasheet

Anti-CSRP2 antibody ab14011

1 References [画像数 1](#)

製品の概要

製品名	Anti-CSRP2 antibody
製品の詳細	Chicken polyclonal to CSRP2
由来種	Chicken
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Mouse, Rat, Human
免疫原	Recombinant fragment, corresponding to amino acids 50-125 of Human CSRP2.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgY

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab14011** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500. Predicted molecular weight: 21 kDa.
IHC-P		Use at an assay dependent concentration. PubMed: 18984963

ターゲット情報

機能	Drastically down-regulated in response to PDGF-BB or cell injury, that promote smooth muscle cell proliferation and dedifferentiation. Seems to play a role in the development of the embryonic vascular system.
組織特異性	Highly expressed in the aorta, but not in heart and skeletal muscle.
配列類似性	Contains 2 LIM zinc-binding domains.
細胞内局在	Nucleus.

画像



E coli-derived fusion protein as test antigen.
 Affinity-Purified IgY dilution: 1/2000, Goat anti-IgY-HRP dilution: 1/1000. Colorimetric method for signal development.

E coli-derived fusion protein as test antigen.
 Affinity-Purified IgY dilution: 1/2000, Goat anti-IgY-HRP dilution: 1/1000. Colorimetric method for signal development.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors