

Product datasheet

Anti-Calretinin antibody [EPR1798(I)] ab76028

リコンビナント RabMAb®

画像数 1

製品の概要

製品名	Anti-Calretinin antibody [EPR1798(I)]
製品の詳細	Rabbit monoclonal [EPR1798(I)] to Calretinin
由来種	Rabbit
アプリケーション	適用あり: WB, IP 適用なし: Flow Cyt, ICC or IHC-P
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide within Human Calretinin (N terminal). The exact sequence is proprietary.
ポジティブ・コントロール	Fetal brain, mouse brain and rat brain lysates.
特記事項	

Our RabMAb® technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMab® patents](#)

This product is a recombinant rabbit monoclonal antibody.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 9% PBS, 40% Glycerol, 0.05% BSA, 50% Tissue culture supernatant
精製度	Tissue culture supernatant
ポリ/モノ	モノクローナル
クローン名	EPR1798(I)
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab76028** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB

IP

追加情報

IP: 1/50.

WB: 1/1000 - 1/5000. Predicted molecular weight: 31 kDa.

Is unsuitable for Flow Cyt, ICC or IHC-P.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報


機能 Calretinin is a calcium-binding protein which is abundant in auditory neurons.

組織特異性 Brain.

配列類似性 Belongs to the calbindin family.

Contains 6 EF-hand domains.

画像



All lanes : Anti-Calretinin antibody [EPR1798(I)] (ab76028) at 1/5000 dilution

Lane 1 : fetal brain lysate

Lane 2 : mouse brain lysate

Lane 3 : rat brain lysate

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : goat anti-rabbit HRP at 1/2000 dilution

Predicted band size: 31 kDa

Observed band size: 31 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors