

## Product datasheet

# Choline Acetyltransferase peptide (736-748) ab45678

### 製品の概要

**製品名** Choline Acetyltransferase peptide (736-748)

### 製品の詳細

**由来** Synthetic

#### アミノ酸配列

**生物種** Human

**配列** CKEKATRPSQGHQP

**領域** 736 to 748

### 特性

Our [Abpromise guarantee](#) covers the use of **ab45678** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**アプリケーション** Blocking

**製品の状態** Liquid

### 前処理および保存

**保存方法および安定性** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### 関連情報

**機能** Catalyzes the reversible synthesis of acetylcholine (ACh) from acetyl CoA and choline at cholinergic synapses.

**関連疾患** Myasthenic syndrome, congenital, 6, presynaptic

**配列類似性** Belongs to the carnitine/choline acetyltransferase family.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors