



## Product datasheet

### Anti-DMRTB1 antibody ab30911

#### 画像数 1

#### 製品の概要

|          |   |
|----------|---|
| 製品名      | Anti-DMRTB1 antibody  |
| 製品の詳細    | Rabbit polyclonal to DMRTB1   |
| 由来種      | Rabbit  |
| アプリケーション | 適用あり: WB  |
| 種交差性     | 交差種: Human  |
| 免疫原      | Synthetic peptide within Human DMRTB1 aa 222-271 (internal sequence). The exact sequence is proprietary.<br>Sequence:<br>LGYLDAPPGVPLQQGFRHVSRSQYQGGGLVSEPGGDFQPSYYLPPPPPL<br><br>Database link: <a href="#">Q96MA1</a>                 |
|          |  <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a> |

#### 製品の特性

|        |   |
|--------|---|
| 製品の状態  | Liquid  |
| 保存方法   | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| バッファー  | Preservative: 0.09% Sodium azide<br>Constituents: 2% Sucrose, PBS                     |
| 精製度    | Protein A purified  |
| ポリ/モノ  | ポリクローナル   |
| アイソタイプ | IgG   |

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab30911** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| アプリ<br>ケーショ<br>ン | Abreviews | 特記事項 |
|------------------|-----------|------|
|------------------|-----------|------|

WB Use a concentration of 2.5 µg/ml. Detects a band of approximately 44 kDa (predicted molecular weight: 36 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

### ターゲット情報

|       |  |
|-------|--|
| 組織特異性 | Testis.  |
| 配列類似性 | Belongs to the DMRT family.<br>Contains 1 DM DNA-binding domain. |
| 細胞内局在 | Nucleus.   |

### 画像



Anti-DMRTB1 antibody (ab30911) at 2.5 µg/ml + HepG2 Cell Lysate at 10 µg

#### Secondary

HRP conjugated anti rabbit IgG at 1/50,000 - 1/100,000

**Predicted band size:** 36 kDa

**Observed band size:** 44 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors