

## Product datasheet

# Anti-Dact2 antibody ab79042

2 References [画像数 2](#)

### 製品の概要

製品名	Anti-Dact2 antibody
製品の詳細	Rabbit polyclonal to Dact2
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB, ELISA, ICC/IF
種交差性	<b>交差種:</b> Mouse, Rat, Human
免疫原	A 12 amino acid synthetic peptide from near the amino terminus of human Dact2 (NP_999627).
ポジティブ・コントロール	SK N SH cell lysate

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.02% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab79042** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 83 kDa.
ELISA		Use at an assay dependent dilution.
ICC/IF		Use a concentration of 20 µg/ml.

## ターゲット情報

機能	Negatively regulates the Nodal signaling pathway, possibly by promoting the lysosomal degradation of Nodal receptors.
配列類似性	Belongs to the dapper family.
ドメイン	The C-terminal PDZ-binding motif may mediate interaction with the PDZ domains of DSH (Dishevelled) family proteins.


## 画像



Anti-Dact2 antibody (ab79042) at 1 µg/ml +  
SK NSH cell lysate at 15 µg

**Predicted band size:** 83 kDa

**Observed band size:** 83 kDa



Immunofluorescence of Dact2 in Human Brain  
cells using ab79042 at 20 µg/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors