

## Product datasheet

# Anti-E. coli O157 antibody (FITC) ab31550

### 1 References

#### 製品の概要

製品名	Anti-E. coli O157 antibody (FITC)
製品の詳細	Goat polyclonal to E. coli O157 (FITC)
由来種	Goat
標識	FITC. Ex: 493nm, Em: 528nm
特異性	Highly specific for E. coli O157:H7
アプリケーション	<b>適用あり:</b> ELISA, ICC/IF
種交差性	<b>交差種:</b> Escherichia coli
免疫原	Whole cells of E. coli serotype O157:H7

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: None Constituents: 50% Glycerol, PBS
精製度	Protein A purified
特記事項(精製)	ab31550 is an affinity purified IgG.
ポリモノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab31550** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		1/20 - 1/2000.
ICC/IF		Use at an assay dependent dilution.

## ターゲット情報

---

### 関連性

*Escherichia coli* are Gram negative bacterium that are commonly found in the lower intestine of warm-blooded organisms (endotherms). Their serological types are determined in combination with somatic antigens (O group: O1-O173) and flagella antigens (H type: H1-H56). The *E. coli* that cause intestinal infectious diseases including diarrhea, acute gastritis or colitis are referred to as pathogenic *E. coli*, which are classified into the following 4 groups according to differences in the mode of pathogenicity; enteropathogenic *E. coli* (EPEC), enteroinvasive *E. coli* (EIEC), enterotoxigenic *E. coli* (ETEC) and enterohemorrhagic *E. coli* (EHEC). Although the identification of pathogenic *E. coli* requires verification of their pathogenicity, pathogenic *E. coli* often have specific serotypes; therefore, typing of the serogroup and serotype is necessary in screening pathogenic *E. coli*.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors