

## Product datasheet

# Anti-ENOSF1 antibody ab102474

1 References [画像数 1](#)

### 製品の概要

製品名	Anti-ENOSF1 antibody
製品の詳細	Rabbit polyclonal to ENOSF1
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	Recombinant protein fragment corresponding to a region within amino acids 63-319 of ENOSF1 (AAH01285).
ポジティブ・コントロール	A549, HeLa, HepG2 or HCT116 cell lysate

### 法規制情報

**医薬用外毒物**

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 20% Glycerol, 1.21% Tris, 0.75% Glycine
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab102474** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 50 kDa.

## ターゲット情報

機能	May regulate thymidylate synthase activity.
配列類似性	Belongs to the mandelate racemase/muconate lactonizing enzyme family, ENOSF1 subfamily.
翻訳後修飾	Could be sumoylated.
細胞内局在	Mitochondrion.

## 画像



Anti-ENOSF1 antibody (ab102474) at 1/1000 dilution + HeLa whole cell lysate at 30 µg

**Predicted band size:** 50 kDa

10% SDS PAGE

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors