


## Product datasheet

# Anti-EXOC6 antibody - N-terminal ab105075

### 画像数 2

#### 製品の概要

製品名	Anti-EXOC6 antibody - N-terminal
製品の詳細	Rabbit polyclonal to EXOC6 - N-terminal
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human <b>交差が予測される動物種:</b> Mouse, Rat, Rabbit, Chicken, Guinea pig, Cow, Cat, Dog, Zebrafish 
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 1-50 (MLEEETDQTY ENVLAEIQSF ELPVEATLRS VYDDQPNAHK KFMEKLDACI) of Human EXOC6 isoform b (NP_001013870). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	HepG2 cell lysate

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab105075** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
----------	-----------	------

WB Use a concentration of 1 µg/ml. Predicted molecular weight: 94 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

### ターゲット情報

**関連性** EXOC6 is a component of the exocyst complex involved in the docking of exocystic vesicles with fusion sites on the plasma membrane.

### 画像




Anti-EXOC6 antibody - N-terminal (ab105075)  
at 1 µg/ml + HepG2 cell lysate at 10 µg

**Predicted band size:** 94 kDa

6-18% gel

Western blot - Anti-EXOC6 antibody - N-terminal (ab105075)



Sample Type: mouse fibroblast lysate (20µg)  
Primary Dilution (ab105075): 1:1000 (2% milk) Secondary Dilution: 1:2000 (5% milk)

Western blot - Anti-EXOC6 antibody - N-terminal (ab105075)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors