

Product datasheet

Anti-Ethinylestradiol antibody ab59670

1 References

製品の概要

製品名	Anti-Ethinylestradiol antibody
製品の詳細	Sheep polyclonal to Ethinylestradiol
由来種	Sheep
アプリケーション	適用あり: ELISA
種交差性	Not applicable
免疫原	Ethinylestradiol conjugated to Bovine thyroglobulin

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.09% Sodium azide Constituents: 0.0536% PBS, 0.87% Sodium chloride
精製度	IgG fraction
特記事項(精製)	Ig fraction prepared by caprylic acid and ammonium sulphate precipitation
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab59670** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use a concentration of 1.25 µg/ml.

ターゲット情報

関連性

Ethinylestradiol, (8R,9S,13S,14S,17R)-17-ethynyl-13-methyl-7,8,9,11,12,14,15,16-octahydro-6H-cyclopenta[a]phenanthrene-3,17-diol, is a semisynthetic alkylated Estradiol with a 17-alpha-ethynyl substitution. It has high estrogenic potency when administered orally, and is often used as the estrogenic component in oral contraceptives. It has a molecular weight 296.40336 g/mol and molecular formula C₂₀H₂₄O₂.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors