

Product datasheet

Anti-FA40B antibody ab103694

画像数 2

製品の概要

|              |   |
|--------------|---|
| 製品名          | Anti-FA40B antibody   |
| 製品の詳細        | Rabbit polyclonal to FA40B  |
| 由来種          | Rabbit  |
| アプリケーション     | 適用あり: WB, IHC-P   |
| 種交差性         | 交差種: Mouse, Human   |
| 免疫原          | Synthetic peptide conjugated to KLH, corresponding to N terminal amino acids 20-50 of Human FA40B |
| ポジティブ・コントロール | Mouse kidney tissue lysates; Human brain tissue   |

製品の特性

|          |   |
|----------|---|
| 製品の状態    | Liquid  |
| 保存方法     | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C long term.             |
| バッファー    | Preservative: 0.09% Sodium Azide<br>Constituents: PBS                                       |
| 精製度      | Immunogen affinity purified   |
| 特記事項(精製) | ab103694 is purified through a protein A column, followed by peptide affinity purification. |
| ポリ/モノ    | ポリクローナル   |
| アイソタイプ   | IgG   |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab103694** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| アプリケーション | Abreviews | 特記事項   |
|----------|-----------|--|
| WB       |           | 1/100 - 1/500. Predicted molecular weight: 86 kDa. |

| アプリケーション | Abreviews | 特記事項   |
|----------|-----------|--|
| IHC-P    |           | 1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. |

## ターゲット情報


**関連性** FA40B belongs to the FAM40 family. Its function is not yet known.

## 画像



Anti-FA40B antibody (ab103694) at 1/100 dilution + Mouse kidney tissue lysates at 35 µg

**Predicted band size:** 86 kDa



ab103694, at a 1/50 dilution, staining FA40B in formalin fixed and paraffin embedded Human brain tissue by Immunohistochemistry, followed by peroxidase-conjugated secondary antibody and DAB staining.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors