


Product datasheet

Anti-FAM62B antibody ab98049

画像数 1

製品の概要

製品名	Anti-FAM62B antibody
製品の詳細	Rabbit polyclonal to FAM62B
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Guinea pig, Cow, Cat, Dog, Pig 
免疫原	Synthetic peptide: corresponding to internal sequence amino acids 611-660 (NSGPNSTIKM KIALRVLHLE KRERPPDHQH SAQVKRPSVS KEGRKTSIKS) of Human FAM62B (NP_065779). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	HT1080 cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab98049** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------


アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 99 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

関連性 FAM62B (Family with sequence similarity 62, member B) is a transmembrane protein which contains three C2 domains. Its specific function is unknown but it may play a role as a calcium-regulated intrinsic membrane protein. There are 5 named isoforms produced by alternative splicing.

細胞内局在 Cell membrane; Multi-pass membrane protein.

画像



Anti-FAM62B antibody (ab98049) at 1 µg/ml
(in 5% skim milk / PBS buffer) + HT1080 cell
lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size: 99 kDa

Gel concentration: 6-18%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors