

## Product datasheet

# Anti-FOXN1 antibody ab113235

1 References 画像数 2

### 製品の概要

製品名	Anti-FOXN1 antibody
製品の詳細	Rabbit polyclonal to FOXN1
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Synthetic peptide, conjugated to KLH, corresponding to residues in Human FOXN1.
ポジティブ・コントロール	Human Uterus tissue; Raji whole cell lysate.

### 法規制情報

医薬用外毒物

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 89% PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab113235** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 69 kDa.
IHC-P		Use a concentration of 10 µg/ml.


## ターゲット情報

機能	Transcriptional regulator involved in development.
組織特異性	Expressed in thymus.
関連疾患	Defects in FOXP1 are the cause of T-cell immunodeficiency congenital alopecia and nail dystrophy (TIDAND) [MIM:601705]. A disorder characterized by the association of congenital alopecia, severe T-cell immunodeficiency, and ridging and pitting of all nails.
配列類似性	Contains 1 fork-head DNA-binding domain.
細胞内局在	Nucleus.

## 画像



ab113235 at 10ug/ml staining FOXN1 in Formalin-fixed, Paraffin-embedded Human Uterus tissue by Immunohistochemistry.



Anti-FOXN1 antibody (ab113235) at 1/1000 dilution + Raji whole cell lysate at 30 µg

**Predicted band size:** 69 kDa

7.5% SDS PAGE.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors