

Product datasheet

Anti-FOXP3 antibody [SP97], prediluted ab99964

リコンビナント RabMAb[®]

4 References 画像数 1

製品の概要

製品名	Anti-FOXP3 antibody [SP97], prediluted
製品の詳細	Rabbit monoclonal [SP97] to FOXP3, prediluted
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human
免疫原	Synthetic peptide within Human FOXP3 aa 400 to the C-terminus (C terminal). The exact sequence is proprietary. Database link: Q9BZS1
ポジティブ・コントロール	Human Tonsil
特記事項	This product is a recombinant rabbit monoclonal antibody.

製品の特性

製品の状態	Prediluted
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.60 Preservative: 0.1% Sodium azide Constituents: 1% BSA, 0.605% Tris
精製度	Protein A/G purified
特記事項(精製)	Purified from TCS by protein A/G.
ポリ/モノ	モノクローナル
クローン名	SP97
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab99964** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		1/1. Boil tissue section in 10mM Citrate, pH 6.0 for 10 min followed by cooling at RT for 20 min.

ターゲット情報

機能	Probable transcription factor. Plays a critical role in the control of immune response.
関連疾患	Defects in FOXP3 are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX) [MIM:304790]; also known as X-linked autoimmunity-immunodeficiency syndrome. IPEX is characterized by neonatal onset insulin-dependent diabetes mellitus, infections, secretory diarrhea, trombocytopenia, anemia and eczema. It is usually lethal in infancy.
配列類似性	Contains 1 C2H2-type zinc finger. Contains 1 fork-head DNA-binding domain.
細胞内局在	Nucleus.

画像



Staining of Human FOXP3 in a Formalin/PFA-fixed paraffin-embedded section of Human Tonsil using prediluted ab99964.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FOXP3 antibody [SP97], prediluted (ab99964)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors