

## Product datasheet

# Goat Anti-Chicken IgY H&L (Alkaline Phosphatase) preadsorbed ab97142

### 1 References

#### 製品の概要

製品名	Goat Anti-Chicken IgY H&L (Alkaline Phosphatase) preadsorbed
由来種	Goat
ターゲット生物種	Chicken
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Chicken IgG and with light chains common to other Chicken immunoglobulins. No antibody was detected against non immunoglobulin serum proteins. Less than 0.1% cross reactivity to bovine, goat, horse, human, mouse, pig, rabbit and rat was detected.
アプリケーション	<b>適用あり:</b> ICC, IHC-P, ELISA, WB
吸着処理血清	Cow, Horse, Human, Mouse, Pig, Rabbit, Rat <a href="#">more details</a>
標識	Alkaline Phosphatase

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, 0.1M Sodium chloride, 50mM HEPES, 1mM Magnesium chloride, 0.1mM Zinc chloride, pH 7.1
精製度	Immunogen affinity purified
特記事項(精製)	Antiserum was cross adsorbed using bovine, horse, human, mouse, pig, rabbit and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Alkaline Phosphatase (AP).
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab97142** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC		Use at an assay dependent dilution.
IHC-P		1/100 - 1/1000.
ELISA		1/2500 - 1/25000. (Primary ELISA)
WB		1/2500 - 1/25000. (Chemiluminescent)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors