

## Product datasheet

# Goat Anti-Chicken IgY H&L (FITC) ab46969

5 References [画像数 1](#)

### 製品の概要

製品名	Goat Anti-Chicken IgY H&L (FITC)
由来種	Goat
ターゲット生物種	Chicken
アプリケーション	<b>適用あり:</b> ICC/IF, IHC-Fr
免疫原	Chicken IgY
標識	FITC. Ex: 493nm, Em: 528nm

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: PBS, 0.2% BSA  Isotonic 0.9%, w/v)
精製度	Affinity purified
特記事項(精製)	Ab46969 is specific to chicken IgY and will not cross react with mammalian IgG.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab46969** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC/IF		
IHC-Fr		

## 追加情報


ICC/IF: 1/250.

IHC-Fr: Use at an assay dependent dilution (PMID 18583367).

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## 画像



[ab4674](#) staining GFAP in human brain cells by Flow Cytometry. Cells were cultured in stem cell media RHB-A and collected using 0.05% Trypsin. The sample was incubated with the primary antibody (1/100 in PBS + 10% FCS) for 20 minutes at 25°C. A FITC-conjugated Goat anti-chicken IgY (ab46969) (1/100) was used as the secondary antibody.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors