

Product datasheet

Goat Anti-Monkey IgG H&L (Alkaline Phosphatase) ab112765

製品の概要

製品名	Goat Anti-Monkey IgG H&L (Alkaline Phosphatase)
製品の詳細	Goat polyclonal Secondary Antibody to Monkey IgG - H&L (AP)
由来種	Goat
ターゲット生物種	Monkey
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Monkey IgG and with light chains common to other Monkey immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. ab112765 may cross react with IgG from other species.
アプリケーション	適用あり: WB, ELISA, IHC-P, ICC/IF
標識	Alkaline Phosphatase

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.10 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, 0.58% Sodium chloride, 0.01% Magnesium chloride, 0.001% Zinc chloride, 1.19% HEPES
精製度	Immunogen affinity purified
特記事項 (精製)	ab112765 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to alkaline phosphatase.
特記事項 (標識)	Molar enzyme: antibody protein ratio is 1:1.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab112765** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/10000. chemiluminescent detection.
ELISA		1/1000 - 1/10000. as a primary.
IHC-P		1/200 - 1/1000.
ICC/IF		1/200 - 1/1000.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors