

Product datasheet

Goat Anti-Mouse IgG F(ab) (TRITC) preadsorbed ab50596

1 References [画像数 1](#)

製品の概要

製品名	Goat Anti-Mouse IgG F(ab) (TRITC) preadsorbed
由来種	Goat
ターゲット生物種	Mouse
特異性	The conjugate can detect mouse IgA, mouse IgM, and all subclasses of mouse IgG. The conjugate shows no reactivity with mouse IgG Fc fragment.
アプリケーション	適用あり: IHC-P, ICC/IF
吸着処理血清	Human more details
免疫原	Purified mouse IgG, Fab fragment
ポジティブ・コントロール	Human peripheral blood lymphocytes and human tonsil sections
標識	TRITC. Ex: 547nm, Em: 572nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 15mM Sodium Azide Constituents: 0.01M PBS, pH 7.4
精製度	Immunogen affinity purified
特記事項 (精製)	The antibody preparation is solid phase adsorbed with human IgG to ensure minimal cross reactivity in tissue or cell preparations.
特記事項 (標識)	Ab50596 is conjugated to crystalline tetramethylrhodamine isothiocyanate (TRITC) in an alkaline reaction, then further purified to remove unbound TRITC.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG


アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab50596** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		1/100.
ICC/IF		1/32.

画像



Immunocytochemistry/ Immunofluorescence - IgG antibody (TRITC) (ab50596)

Immunofluorescent staining of microglia in the hippocampus (CA3 region) in rat brain. Rat brain sections stained with Monoclonal Anti-Glial Fibrillary Acidic Protein (clone no. G-A-5), and ab50596 at 1/32 dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors