

Product datasheet

Goat Anti-Mouse IgG H&L preadsorbed ab98800

製品の概要

製品名	Goat Anti-Mouse IgG H&L preadsorbed
製品の詳細	Goat polyclonal Secondary Antibody to Mouse IgG - H&L, pre-adsorbed
由来種	Goat
ターゲット生物種	Mouse
アプリケーション	適用あり: IHC-P, ELISA, WB, ICC
吸着処理血清	Human, Rat

[more details](#)

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
特記事項(精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab98800** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		1/200 - 1/2000.
ELISA		1/100 - 1/30000. Coating: Use 1/100-1/500 dilution; Primary: Use 1/1,000-1/30,000 dilution.

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/30000. Colorimetric: Use 1/1,000-1/10,000 dilution; Chemiluminescent: Use 1/1,000-1/30,000 dilution.
ICC		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors