

Product datasheet

Anti-HYAL2 antibody ab90004

画像数 2

製品の概要

製品名	Anti-HYAL2 antibody
製品の詳細	Rabbit polyclonal to HYAL2
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	A peptide near intermediate residues of Human HYAL2 protein (NM_003773).
ポジティブ・コントロール	Human fetal thymus; SW480 cell lysate

法規制情報

医薬用外毒物

製品の特性

製品の状態	Lyophilised:Reconstitute with 200ul distilled sterile water. Please note that if you receive this product in liquid form it has already been reconstituted as described and no further reconstitution is necessary.
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 2% BSA, pH 7.2
精製度	Immunogen affinity purified
特記事項(精製)	Purified by peptide affinity column
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab90004** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/200 - 1/1000. Predicted molecular weight: 54 kDa.
IHC-P		1/50 - 1/200.


ターゲット情報

機能	Hydrolyzes high molecular weight hyaluronic acid to produce an intermediate-sized product which is further hydrolyzed by sperm hyaluronidase to give small oligosaccharides. Displays very low levels of activity. Associates with and negatively regulates MST1R.
組織特異性	Widely expressed. No expression detected in adult brain.
配列類似性	Belongs to the glycosyl hydrolase 56 family. Contains 1 EGF-like domain.
細胞内局在	Cell membrane.

画像



ab90004, at a dilution of 1/100, staining HYAL2 in the membrane of formalin fixed, paraffin embedded Human fetal thymus by Immunohistochemistry.
(A): No immunizing peptide
(B): Blocked with immunizing peptide



Anti-HYAL2 antibody (ab90004) at 1/500 dilution + SW480 cell lysate

Predicted band size: 54 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors