

## Product datasheet

# Anti-Helicobacter pylori Polysaccharide Deacetylase antibody ab93827

### 製品の概要

製品名	Anti-Helicobacter pylori Polysaccharide Deacetylase antibody
製品の詳細	Rabbit polyclonal to Helicobacter pylori Polysaccharide Deacetylase
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Recombinant fragment 交差が予測される動物種: Other species
免疫原	Synthetic peptide corresponding to N terminal residues of <i>Helicobacter pylori</i> P12 Polysaccharide Deacetylase (UniProt ID: B6JKN9).

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	Preservative: 0.01% Sodium Azide Constituents: 50% Glycerol, PBS
精製度	Immunogen affinity purified
特記事項(精製)	purity >90%
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab93827** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項

アプリ ケー ション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 34 kDa. for 2 hours. This antibody has been tested in Western blot against the recombinant peptide used as an immunogen. We have no data on detection of endogenous protein.
ELISA		1/2000 - 1/5000.

## ターゲット情報

<b>関連性</b>	Polysaccharide Deacetylase contains CDA1 which is predicted xylanase/chitin deacetylase [Carbohydrate transport and metabolism] and domain of Polysacc deac 1; Polysaccharide deacetylase. Polysaccharide Deacetylase functions in biosynthesis of surface polysaccharides and lipopolysaccharides. The function proposed based on presence of conserved amino acid motif, structural feature or limited homology; putative enzyme.
------------	---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors