

Product datasheet

Anti-Hemopexin antibody [4] (Biotin) ab121081

製品の概要

製品名	Anti-Hemopexin antibody [4] (Biotin)
製品の詳細	Mouse monoclonal [4] to Hemopexin (Biotin)
由来種	Mouse
標識	Biotin
アプリケーション	適用あり: Sandwich ELISA
種交差性	交差種: Human
免疫原	Hemopexin purified from Human plasma adsorbed onto aluminum hydroxide gel.
特記事項	Myeloma: SP2mIL6

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40 Preservative: 0.1% Sodium azide Constituents: 99% Phosphate Buffer, 0.82% Sodium chloride
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	4
アイソタイプ	IgG2b
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab121081** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Sandwich ELISA		Use at an assay dependent concentration. Can be paired for Sandwich ELISA with Mouse monoclonal [32] to Hemopexin (ab27711) . ab121081 works as a detecton antibody.

ターゲット情報

機能	Binds heme and transports it to the liver for breakdown and iron recovery, after which the free hemopexin returns to the circulation.
組織特異性	Expressed by the liver and secreted in plasma.
配列類似性	Belongs to the hemopexin family. Contains 5 hemopexin-like domains.
翻訳後修飾	N- and O-glycosylated. O-glycosylated with core 1 or possibly core 8 glycans.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors