

## Product datasheet

### Hsp90 peptide ab40009

#### 製品の概要

---

製品名 Hsp90 peptide

#### 法規制情報

---

医薬用外毒物

#### 製品の詳細

---

由来 Synthetic

#### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab40009** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Western blot  
Blocking - Blocking peptide for Anti-Hsp90 antibody ([ab19021](#))

製品の状態 Liquid

#### 前処理および保存

---

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.  
Preservative: 0.02% Thimerosal (merthiolate)  
Constituents: 0.1% BSA, PBS, pH 7.2

#### 関連情報

---

機能 Molecular chaperone that promotes the maturation, structural maintenance and proper regulation of specific target proteins involved for instance in cell cycle control and signal transduction. Undergoes a functional cycle that is linked to its ATPase activity. This cycle probably induces conformational changes in the client proteins, thereby causing their activation. Interacts dynamically with various co-chaperones that modulate its substrate recognition, ATPase cycle

	and chaperone function.
<b>配列類似性</b>	Belongs to the heat shock protein 90 family.
<b>ドメイン</b>	The TPR repeat-binding motif mediates interaction with TPR repeat-containing proteins like the co-chaperone STUB1.
<b>翻訳後修飾</b>	ISGylated. S-nitrosylated; negatively regulates the ATPase activity and the activation of eNOS by HSP90AA1.
<b>細胞内局在</b>	Cytoplasm. Melanosome. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors